



THE EXTREMITIES SOLUTION




Orfit offers a complete immobilization solution, specifically designed for fixation of extremities in a reproducible way.


The system consists of a low density carbon fiber base plate with a leg support, two comfort cushions and dedicated masks. The base plate can be indexed to any couch top on 2-pin bars and has multiple holes to offer a maximum of possibilities for positioning and immobilizing both upper and lower extremities.

When using the leg support, one leg can be lifted and placed out of the treatment field. Masks can be made on the leg support as well as on the base plate itself. Foam cushions are placed under the knee and heel to provide patient comfort. Upper extremities can be immobilized while the patient is lying on the couch, or sitting next to it. In both cases, masks are made over the elbow and hand to obtain an immobilization with excellent reproducibility. Comfort cushions are used under the elbow and the hand for patient comfort.






EXTREMITIES BASE PLATES




Art. N° 32103
Set - CFL



Art. N° 32103/1
Base plate - CFL




Art. N° 32103/3
Leg support - CFL




Art. N° 32103/MR
Set - MR compatible


ANTIBACTERIAL MASKS FOR EXTREMITIES




Art. N° 35792/32MA
3-points foot mask



Art. N° 35793/32MA
2-points knee mask



Art. N° 35794/2MI
2-points hand mask



Art. N° 35795/2MI
2-points elbow mask - straight



Art. N° 35796/2MI
4-points elbow mask - curved



Art. N° 35799/2MI/R
5-points shoulder mask (right shoulder)

Art. N° 35799/2MI/L
5-points shoulder mask (left shoulder)